INSTRUCTIONS FOR USE

CyberCore

Radioopaque Composite core build-up material

- 1.) Isolate the tooth from oral fluids. A rubber dam is recommended.
- 2.) Acid-etch the entire cavo surface with a 32% to 40% phosphoric acid (like for instance CyberEtch) for 15 seconds. (If tooth is to be prepared for a post, etching of canal and cavo surface can be completed simultaneously after post is sized and fitted for canal).
- 3.) Rinse off etchant with a water-spray for at least 15 seconds. Etched enamel should appear frosted.
- 4.) Blot-dry etched surfaces.
- 5.) Follow bonding agent instructions for use with dual-cure composites.
- 6.) Select shade and assemble syringe mixtip. (Bleed syringe on initial use and bleed through assembled mixing tip and intraoral tip before each use).
- 7.) Place intraoral tip of assembled syringe into base of prepared tooth and dispense an adequate amount of material to complete build-up. Let material rise around the tip while withdrawing to eliminate air bubbles.
- 8.) Do not disturb during the gelation stage. Light-cure for at least 40 seconds. Check for complete hardness. Chemical-cure intra-oral set-time is approximately 4 minutes from start to mix/application. Set-time takes longer outside of mouth.
- 9.) Finish preparation with a bur.

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