**GB - CyberSilicone Vinylpolysiloxane (Hydro) Impression Material. Light Body, Fast Set and Light Body, Regular Set, Regular Body INSTRUCTIONS FOR USE**

**Description:**  
ISO 4823 Type 3 Light-bodied Consistency - ADA 19 Type 1 Low Viscosity. ISO 4823 Type 2 Medium-bodied Consistency - ADA 19 Type 1 Medium Viscosity.

CyberSilicone Vinylpolysiloxane (Hydro) is a hydrophilic addition reaction silicone impression material dispensed in an automix cartridge delivery system. CyberSilicone Hydro Impression Material is thixotropic - it resists slumping, yet flows under pressure. CyberSilicone Hydro Light Body is available in two setting types: Fast Set and Regular Set. Light Body, Fast Set is designed primarily for use with 1 or 2 step impressions where a faster setting time is required.

**Indications:** Correction impression material for the two-step (putty-wash) impression technique. Use with dispenser to extrude. Use Intra-oral Tips for extruding material directly into mouth.

**Impression Procedure:** Inject automix material evenly around prepared teeth and/or directly onto Putty Impression. Seat tray and maintain in position for approximately 2,30 minutes (Fast Set) or 3,30 minutes (Regular Set). **Note:** Set time is based on 23°C. Higher temperature decreases set time, lower temperature increases set time.

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**Electroplating:** Impressions may be copper or silverplated.

**Storage:** Do not store above 27°C (80°F). Shelf life 3 years when stored at 23°C (73°F).

**General Notes:** Do not use with condensation reaction silicone impression materials.

TECHNICAL DATA	REGULAR SET	FAST SET
Mixing Ratio Base : Catalyst	1 : 1	1 : 1
Mixing Time	Auto Mixing	Auto Mixing
Total Working Time* (incl. Mixing Time)	2 min.	1 min. 30 secs.
Time in mouth* (minimum)	3 min. 30 secs.	2 min. 30 secs.
Setting Time* (includes Mixing time)	5 min. 30 secs.	4 min.
Linear Dimensional Change after 24 Hrs.	<0,20%	<0,20%
Maximum Strain in Compression	3%-5%	3%-5%
Recovery from Deformation	>99,5%	>99,5%
Hardness « Shore A »	45±2	45±2

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Light Body, Fast Set  
Light Body, Regular Set  
Regular Body

Hydrophilic/Hidrófilo/Hydrophile/Hydrophil/Hidrofilo

Vinylpolysiloxane Impression Material, Light Body, Fast Set, Light Body Regular Set & Regular Body  
• INSTRUCTIONS FOR USE

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DE Healthcare Products  
Gillingham  
MEB OSB U.K.

## GB - CyberSilicone Vinylpolysiloxane Putty INSTRUCTIONS FOR USE

**Description:**

ISO 4823 Type 0 Putty Consistency - ADA 19 Type 1 Very High Viscosity. CyberSilicone Vinylpolysiloxane Hydro putty is a hydrophilic addition reaction silicone impression material.

**Indicaciones:** Preliminary impression for the two-step (putty-wash) technique in combination with CyberSilicone Hydro Light Body, Regular Set/Fast Set or Regular Body syringe materials.

**Trays and Adhesives:** Perforated or non-perforated metal or plastic trays may be used with CyberSilicone Hydro putty. To improve adhesion of Putty

to the tray, adhesive should be applied in a thin, even layer. Allow adhesive to dry for 3 minutes before loading the impression tray.

**Dispensing and Mixing:** Dispense equal volumes or weights of catalyst and base, using the colored spoon for the base and the white spoon for the catalyst. Do not interchange spoons. Knead with fingertips until uniform and streak-free (approximately 30 seconds). **Note:** Containers should be closed immediately after use, taking care not to interchange lids.

**Impression Procedure:** Load properly mixed material into the prepared tray. Place a polyethylene spacer sheet over putty and seat tray in the mouth not longer than 1 minutes after start of mix. Maintain position for approximately 2 minutes (material should have a rubbery consistency). **Note:** Higher temperature decreases set time; conversely, lower temperature increases set time.

**Casting:** Store impression at room temperature. Impression may be poured 30 minutes after removal from mouth.

**Electroplating:** Impressions may be copper or silver plated.

**Storage:** Do not store above 27°C (80°F). Shelf life is 3 years when stored at 23°C (73° F).

**General:**

- Do not interchange catalyst and base covers. Do not interchange measuring spoons.
- Certain gloves will inhibit the setting of putty. Test gloves to confirm proper setting times. Vinyl gloves or vinyl overliners are recommended. DO NOT USE LATEX GLOVES.
- Do not use with condensation reaction silicone impression materials.

TECHNICAL DATA	REGULAR SET	FAST SET
Mixing Ratio Base <span> </span> : Catalyst	1 <span> </span> : 1	1 <span> </span> : 1
Mixing Time	30 secs.	30 secs.
Total Working Time* (incl. Mixing Time)	2 min.	1 min. 30 secs.
Time in mouth* (minimum)	3 min. 30 secs.	2 min. 30 secs.
Setting Time* (includes Mixing time)	5 min. 30 secs.	4 min.
Linear Dimensional Change after 24 Hrs.	<0.20%	<0.20%
Maximum Strain in Compression	1%-3%	1%-3%
Recovery from Deformation	>99%	>99%
Hardness « Shore A »	60±2	60±2

\*23°C/73°F

**Disinfection:** Impression may be disinfected or sterilized with a disinfectant for impressions.

**Compatibility with Gypsum:** This product is compatible with all well known gypsum materials.

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