CE 0473

DE Healthcare Products
Gillingham
ME8 0SB U.K.

CYBER TECH

FILL

NANO HYBRID
COMPOSITE

Nano-Hybrid Composite • Resin-based dental restorative material

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ISO 4049, Type 1, class 2, group 1, radio-opaque

DE Healthcare Products
Gillingham
ME8 0SB U.K.

INSTRUCTIONS FOR USE

CyberFill Nano Hybrid Composite

Resin based dental restorative material

Cybertech Fill Nano-Hybride Composite is a syringeable, sculptable, radiopaque, highly filled nano-hybrid composite for the restoration of all teeth. Cybertech Fill Nano-Hybride Composite covers a wide range of tooth shades with few shades.

Technical Data

• Complies with ISO 4049:2000 • Average filler particle size: 0.4 µm

Cybertech Fill Nano-Hybrid

• Percentage by volume of total inorganic filler: 57% • Percentage by weight of total inorganic filler: 76%

Indications

Cybertech Fill Nano-Hybride Composite is specifically indicated for:

• Direct filling of all teeth • More translucent material in the layer technique • Bonding of composite and ceramic restorations

Contra-indications

If the site cannot be isolated after enamel etching and during application and curing of Cybertech Fill Nano-Hybride Composite. If allergic to components of Cybertech Fill Nano-Hybride Composite. If oral hygiene is poor.

Side effects

No systemic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.

Interactions with other agents

Agents which contain eugenol and/or oil of cloves may affect the polymerization of Cybertech Fill Nano-Hybride Composite. The use of zinc-oxide-eugenol cements in combination with Cybertech Fill Nano-Hybride Composite should therefore be avoided. Discolorations may occur when using cationic mouth rinses, plaque indicators or chlorhexidine.

Application - syringe

Cleansing the tooth. The tooth to be treated and the adjacent teeth are cleansed using a brush and prophylaxis-paste which does not contain fluoride. Selection of shade. Select shades before isolation. Ideally, the shade is determined in broad daylight using the VITA™ Lumin-Vacuum shade guide (register-red trademark of the company VITA™). Shades are identical to those of the shade guide only 24 hours after curing. Isolation of the cavity (rubber dam) Adequate isolation is absolutely necessary for achieving optimum results. Preparation of the cavity:

When preparing the cavity, every effort should be made to use a technique that is gentle to the tooth substance (principle of Adhesive Restoration Technique). Prepare the enamel and dentine using 80 µm preparation diamonds and finish using 25 µm finishing diamonds. Bevelled enamel margins are recommended to increase the area of adhesion between tooth and filling material, thereby optimizing marginal seal. Pulp Protection: Covering the areas close to the pulp with hardsetting calcium hydroxide cements using the spot technique will provide protection of the pulp against bacterial infiltration. Bonding system: use a bonding according to the corresponding instructions for use.

Matrix and interdental wedges Apply a thin matrix band for class II, III and IV restorations. Fix the matrix proximally using interdental wedges. **Application of Cybertech Fill Nano-Hybride Composite from a syringe:** Dose Cybertech Fill Nano-Hybride Composite directly onto the mixing pad, take up in small increments with a composite instruments or any other application instrument, and apply to the cavity.

Application – tip: Remove the cap. Insert the composite capsules through the barrel opening of the dispenser. Rotate the angled composite capsule into the desired position. Extrude Cybertech Fill Nano-Hybride Composite into the cavity or onto an instrument by exerting gentle, even pressure on the handle. Working time Cybertech Fill Nano-Hybride Composite is photosensitive, and should not be exposed to light too long, in particular to the operating light or sunlight. If possible, protect exposed Cybertech Fill Nano-Hybride Composite material by means of light opaque covering (not blue).

Polymerization

Cybertech Fill Nano-Hybride Composite can be cured by using current dental light sources. Recommended exposure times: 40 sec. per layer (2mm). When multiple surfaces are available, cure each surface 40 sec. Do not remove the inhibition layer since this would affect adhesion between the layers. In the event of underexposure there is a risk that the ideal strength values may not be achieved. Overexposure is not possible. Equipment whose light output is unknown should be checked with a light meter according to corresponding instructions for use.

Finishing procedures

Finishing consists of three stages: coarse finishing, fine finishing and polishing. Finishing diamonds (40 µm and 15 µm) are the most versatile and the least destructive. They should be used with very light pressure, with a constant wiping motion and a generous water spray. The ideal speed ranges between 5000 and 15 000 rpm. Flexible diamond files (80 µm, 40 µm and 15 µm) and strips are suitable for aproximal surfaces. For high gloss polishing results use any comparable instruments. Fluoridize all treated surfaces after finishing.

Emergency measures

Danger from swallowing of Cybertech Fill Nano-Hybride Composite can be excluded due to the small quantity applied and its lack of toxicity. In case of direct contact with the oral mucosa, rinsing with tap water is sufficient. In the case of contact with the eyes, rinse thoroughly with water (10 min) and consult physician.

Toxicological characteristics

In addition to the active components, Cybertech Fill Nano-Hybride Composite contains small amounts of stabilizers, activators and color pigments. Allergic reaction to resin components or the amine component of the activator system is possible.

Note

Only supplied to dentists and dental laboratories or upon their instructions. Keep out of the reach of children! Do not mix different shades on the mixing pad, since this procedure will lead to air entrapment and porosities.

Storage conditions and shelf life:

Store at 5–27 °C / 41–80 °F. The expiration date is marked on the packaging. Do not use after expiration date. Every effort should be made to protect the composite from visible light prior to use.

Attention

Composite capsules are designed for single patient use only. Do not re-cap and / or re-use the capsules once material has been dispensed for that patient.

Caution

Federal law restricts this device to sale by or on the order of a dentist. Uncured methacrylate resin may cause dermatitis and damage the pulp. Avoid contact to skin, eyes and soft tissues. Wash thorough with water after contact.

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