CYBER TECH

FILL 45

NON GAMMA II HIGH COPPER

Amalgam Non Gamma 2. High copper 45%.

Dental filling material.

• **INSTRUCTIONS FOR USE**

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Alloy for producing non-gamma 2 amalgam.

Active ingredient: non-gamma 2 amalgams

Composition of the medical product: Dental filling material

Active ingredients: amalgams are produced by mixing approximately equal weights of alloy powder and mercury to create a plastic material which the dentist applies. The powders consist of lathe-cut and spherical particles.

Composition of the alloy powder – 100 g contain:

Ag: 45.5% / Sn: 31.5% / Cu: 23%

After hardening, amalgams mostly consist of Ag, Sn, Cu3Sn original phases, frequently Cu/Ag eutectic alloy, Ag3Hg4 and Cu6Sn5 reaction products and no, or only a minimal amount, of Sn8Hg. As the efficacy of dental materials depends heavily on the pharmaceutical/material quality, these instructions were compiled for non-gamma 2 amalgams – their properties comply with ISO 24234:2004 Mercury and alloys for dental amalgam.

Range of indications: Amalgam must only be used for occlusal load-bearing, posterior fillings (classes I and II) and only if other composite filling materials cannot be used. As a preventive measure, the number of amalgam fillings per patient must be kept to a minimum as each amalgam filling exposes the body to mercury. **Contraindications:** Contraindications are illnesses or conditions where certain medical products must not be used at all or only after careful consideration by the dentist as, in general, the anticipated effect relates unfavourably to the possible harm. To enable the dentist to check carefully for contraindications, he must be informed about previous and concomitant illnesses, any simultaneous treatment and any special circumstances and habits. Contraindications may not occur or become known until after commencing treatment with this medical product. You should also inform your dentist in these cases. Patients with severe functional disorders of the kidneys should not be provided with new amalgam fillings. Hypersensitivity (allergy) to amalgam constituents. Amalgam is not indicated: for

retrograde root fillings, for building up cores for supporting crowns and inlays or for filling cast crowns – new amalgam fillings must not be placed in occlusal or approximal contact with existing cast restorations. Due to their increased sensitivity to mercury, all children below the age of 6 and particularly those younger than 3, should only be treated with amalgam after very careful consideration. **Side-effects:** Apart from the desired effects, pharmaceuticals may also cause undesirable effects known as side-effects. The following are side-effects observed while CyberFill was in use, but which may not necessarily affect every patient. In specific cases, hypersensitivity (type IV allergic reactions) and/or electrochemical reactions (e.g. taste sensations and irritation, i.e. of the oral tissues) have been described. After an amalgam filling has been placed or removed, the mercury concentration in the blood and urine increases temporarily. Electrochemical processes may provoke the lichen planus (eczema) of the oral mucosa.

Interactions with other materials: The effects of some medical products may be influenced by using other materials simultaneously. Therefore, if you take other medication continually, have taken it until just recently or wish to take it simultaneously with this medical product you should ask your dentist. Your dentist will inform you whether incompatibilities have to be reckoned with under these circumstances or if special measures such as a change in dosage are required for you to use this medical product. If the occlusal or other surfaces contact opposing crowns, bridges, inlays or fillings made of other alloys, galvanic effects (caused by corrosion) may occur. If persistent local malaises (e.g. changes in taste sensation) occur with amalgam fillings in contact with other alloys due to electrochemical effects, the amalgam filling must be replaced with other materials. This applies to the following clinical situations: core material for supporting fixed restorations, filling material in crowns made of other alloys, filling material in approximal or occlusal contact with other alloys.

Warnings: Placing an amalgam filling without protecting the pulp/dentine adequately, especially in deep or medium cavities which have not been lined, may irritate the pulp. Condensing amalgam into the gingival sulcus of compound cavities injures/destroys the periodontium in this region. Should lichen planus occur on oral tissues in contact with amalgam, it is advisable to replace these amalgam fillings with non-metallic materials. Lichen planus is a polyaeitiological dermatosis and in its basic form does not require treatment in the oral cavity. With specific types of lichen planus, replacing metals may be successful. Lichen planus may be provoked by an isomorphous stimulus, which may include electrochemical processes.

Important incompatibilities: None known. **Dispensing:** The following guidelines apply unless otherwise prescribed: the amount of amalgam mixed depends on the size of the cavity. The dentist determines the amount. **Type and duration of use:** CyberFill is mixed (trituated) with an ordinary amalgamator: **CyberFill capsules:** 4 - 5 secs. Once the cavity has been isolated (e.g. with a rubber dam), portions of the soft material are placed in it and condensed with adequate pressure. When condensing the approximal regions of the cavity, it is essential to use a suitable matrix. The cavity must be slightly overfilled at first before removing the mercury rich excess amalgam and contouring the filling. **Working time:** Condense 3.5-5 min, Band removal, carving and burnishing 5-9 min, Polish 24 hr after placement. **Setting time:** 3.5 - 5 min. The filling must be polished one day later at the very earliest. Cool by spraying continually and adequately. Check the occlusion with articulating tool. Taking the correct precautions when removing amalgam fillings reduces the exposure of the patients and staff – these precautions include using an aspirator, placing a rubberdam, spraying adequate water, ventilating the surgery and disposing of waste amalgam according to regulatory requirements. The dentist decides on the duration of treatment. **Overdosage:** Not applicable.

First aid, symptoms and antidote: None. **Pharmacological and toxicological properties, pharmacokinetics, bio-availability Pharmacological properties:** After hardening, amalgam is an alloy of mercury and other metals consisting of various intermetallic phases. Amalgam fillings release elementary mercury and constituents of Cu6Sn5, Sn8Hg, Sn(OH)Cl and CuCl2 x 3 Cu(OH)2. **Acute toxicity:** While placing amalgam fillings, during the first few days thereafter and when removing amalgam fillings increased amounts of mercury vapour are released. Increased concentrations of mercury may be recorded temporarily in the saliva, blood and urine. **Chronic toxicity:** In case of chronic uptake of inorganic mercury compounds, the brain, liver and especially the kidneys are considered critical organs. The mercury concentrations found in these organs and the mercury content in the saliva, blood and urine correlate with the number of amalgam fillings, respectively, filling surfaces. Various animal experiments proved the presence of mercury in these organs after placing amalgam fillings. High mercury values (e.g. 1.86 mg/kidney in sheep) did not correlate with the clinical symptoms. Examinations of the kidney tissue of corpses indicated that persons (n = 7) with 11 – 33 amalgam filling surfaces have a higher mean mercury concentration (433(48 – 810)ng Hg/g (wet weight)) than those without amalgam filling surfaces (49(21 – 105) ng Hg/g (wet weight)). Examinations of corpse brains showed that the mercury concentration correlates with the number of amalgam filling surfaces: Persons (n = 51) with 5 – 14.5 amalgam filling surfaces exhibited a mean mercury concentration of 1521 (3 – 121.4 ng/g brain (grey substance)) whereas persons with 0 – 1 amalgam filling surfaces exhibited a mean mercury concentration of 6.7 (1.9 – 22.1) ng/g brain. Persons with amalgam fillings were shown to have higher mercury concentrations in their blood and urine than persons without amalgam fillings: Where more than 10 amalgam fillings were present, the geometrical mean value was 0.61 lg Hg/l blood, respec-

tively, 1.45 lg Hg/ urine. In comparison, persons with no amalgam fillings exhibited mercury concentrations of 0.46 lg Hg/l blood, respectively, 0.28 lg Hg/l urine. The WHO (1976) quotes a lowest value of 35 lg Hg/l blood, respectively, 150 lg Hg/l urine above which one can reckon with a mercury effect in sensitive persons with unspecific symptoms. As a preventive health measure, the WHO suggests (1980) a value of 50 lg Hg/l urine. The first biological effects in the kidneys (such as increased enzymuria and proteinuria) are observed above 50 lg Hg/l urine. In 1991, the WHO determined that exposure to mercury vapour within a range of 25 – 80 lg/m3, which corresponds to a urinal mercury content of 30 – 100 lg Hg/l creatinine, "increases the incidence of certain less serious toxic effects which do not lead to manifest clinical impairment". Thus, the lowest observed effect level (LOEL) of mercury concentration in urine is 30 lg Hg/g creatinine. According to the evaluation criteria defined by the (German) Federal Health Authorities in 1987, health hazards cannot be ruled out above a mercury concentration of 10 lg Hg/l blood, respectively, 20 lg Hg/l urine (KRAUSE, C. et. al. in: "Schr.-Reihe Verein WaBoLu 74, 105 (1987)).

Reproduction toxicology: Mercury and its compounds are transferred via the placenta and small amounts enter the mother's milk. Post mortems on humans showed that the level of mercury in the organs of fetuses (e.g. liver) and new-born babies (e.g. kidneys) is linked to the number of amalgam fillings in the mother's mouth. However, no macroscopic changes in the tissue of these organs were found. Animal experiments and epidemiological tests also produced no founded suspicion that this exposure to mercury leads to prenatal or postnatal health hazards. **Pharmacokinetics:** The following information refers to mercury released in the form of metallic (liquid or vaporized) and ionized mercury while amalgam fillings were being processed, were in situ or being removed. The causes of interindividual differences in mercury release and resorption in principally the same amalgam restoration remain to be investigated. The WHO estimates the amount of mercury released from amalgam fillings to be 3.8 lg to 21 lg per day – between 3 lg and 17 lg of which is resorbed per day. Thus, the additional amount of mercury taken up from amalgam fillings is calculated as being a maximum of 6.5 times the average uptake of 2.61 lg per day from other sources (WHO (1991)). The uptake of metallic mercury in the digestive tract is less than 1%. The gastrointestinal resorption rate for mercury ions (Hg22+ and Hg2+) is between 7% and 10%. Approximately 80% of vaporized mercury is resorbed via the lungs. This is no adequate proof that mercury vapour is transported directly via the nasal mucosa to the brain. Mercury is oxidized, partially reduced and possibly methylated in small amounts by the organism. Mercury has an affinity for compounds such as albumin, glutathione or cystein which support sulphydryl groups. Mercury is linked to metallothionein metallic bonding protein in various organs and accumulates. Ionized mercury is not soluble in lipids and will not penetrate the blood-brain barrier. It concentrates itself in parenchymatic organs such as the liver and especially the kidneys. Mercury absorbed in vapour form is lipophilic and passes the blood-brain barrier. It may accumulate in parts of the brain. Ionized mercury or that absorbed in vapour form passes the placenta barrier. The half-life of mercury in the organism and kidneys is approximately 60 days, but can exceed one year in the brain. Mercury is excreted renally, sterocorally and with exhaled air. Cleaning the teeth and chewing gum increases the mercury values in exhaled air and saliva significantly. As different test methods were used, it is difficult to compare the values. There are indications that so-called amalgam tattoos, caused by amalgam particles or corrosion products, on the gingiva or oral mucosa may release mercury depending on the size of the particles.

Pregnancy and lactation period: As the foetus is exposed to mercury from the mother's amalgam fillings, as a preventive measure no, or no more, amalgam should be placed during the pregnancy. Wherever possible, alternative materials should be given preference. As additional mercury is released by removing amalgam fillings, clinically impeccable amalgam fillings should not be removed, especially during pregnancy. Current knowledge provides no proof that exposing the foetus to mercury from the mother's amalgam fillings has an effect on the child's health. **Effects on drivers and operating machinery:** None. **Shelf-life:** The shelf-life is 5 years. Do not use after the expiry date. **STORAGE** **For dentists:** Hazards to dentists and staff caused by Hg vapour can be avoided by handling the Hg and amalgam correctly during processing, keeping waste amalgam in fixing salt solution in tightly sealed containers and ventilating the practice well. **Special precautions for disposing of unused amalgam:** Comply with the relevant waste disposal regulations. **Presentation and packaging** **CyberFill capsules, 1 spill:** 400 mg alloy capsules **CyberFill capsules, 2 spill:** 600 mg alloy capsules **CyberFill capsules, 3 spill:** 800 mg alloy capsules **Keep out of reach of children!** **For use by dental professionals only!**

1) These geometrical mean values correspond to the as yet unpublished results of the 2nd environmental survey on mercury for the Federal Republic of Germany, 1990/92 (Institute for Water, Ground and Air Hygiene – Federal Ministry of Health).

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