# **INSTRUCTIONS FOR USE**

## **CyberCement Resin**

**Dual cured permanent adhesive resin cement**

### **Ingredients:**

#### **Base:**

- Aliphatic polyester urethane acrylate uretano
- Hydroxyethyl methacrylate • Crosslinking monomer
- Photoinitiators • Adhesive accelerator • Glass filler

#### **Catalyst:**

- Triethyleneglycol dimethacrylate
- Bis-phenol-A glycidildimethacrylate
- Benzol peroxide • Glass filler

### **Indications:**

CyberResin Cement is a multi-purpose, dual cured, radiopaque, permanent cement of crowns and bridges, inlays and onlays, posts and cores, ceramic crowns and Maryland bridges, Veneering of alloys, Implant prosthesis. It can be used with base alloys, porcelain, dentine, enamel and with amalgams. It is a self cured cement which can also be used with light. Light curing of CyberResin Cement allows a better strength at the margins more stability.

### **Setting time:**

Working time at room temperature is at least 2 minutes.

Final setting time at room temperature is approximately 5 minutes.

Oral set time is approximately 3,5 minutes. Light cure for 20 seconds.

### **Caution:**

CyberResin Cement contains polymerisable monomers which may cause allergic reactions with susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs or if known allergy to methacrylate resin exists, discontinue use.

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