INSTRUCTIONS FOR USE

CyberSeal

Light cured, Pit & fissure sealant with fluoride release.

- *For the preventive sealing of pits and fissures.*
- *Easy & precise to handle, low viscosity, excellent flow properties, optimal wetting of fissures.*
- *Fluoride release.*

INSTRUCTIONS FOR USE:

- Clean the enamel surface thoroughly with either cleansing paste containing neither oil nor fluoride, or with pumice and brush. When preparing fractured porcelain crowns for bonding, abrade the surface first.
- Isolate the working field, preferably with a rubber dam
- Apply the CyberEtch gel, and let it work for 30 to 60 seconds
- Rinse thoroughly
- Dry with air free of water and oil. The etched enamel should have a mat white appearance. Avoid contamination of the etched surface with saliva
- Mount and secure lock a new dispensing tip on the CyberSeal syringe and carefully syringe an adequate amount of it into the fissures. Let the sealant flow until the entire fissure is covered. You can use a disposable brush or the dispensing tip itself to spread the sealant and eliminate air bubbles.
- Cure sealant for 20 sec
- Check occlusion

INGREDIENTS:

- Urethandimethacrylate oligomer
- Bis-GMA
- Triethyleneglycoldimethacrylate
- 2-Hydroxyethylmethacrylate
- Glass fillers
- Photoinitiators

CONTRAINDICATIONS:

- If pits or fissures are already affected with caries
- If a dry working field cannot be established during the application and curing of CyberSeal.
- If the patient is known to be allergic to any of the ingredients of CyberSeal CyberSeal contains polymerizable monomers, which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

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