INSTRUCTIONS FOR USE

CyberTemp Long

CyberTemp Long – Long term temporary & implant cement, non eugenol. The material can be removed easily. It is stronger than a regular temporary cement.

- Check occlusion and proximal contacts. Complete all finishing procedures.
- Dry and isolate the prepared abutments. Absolute dryness is not necessary, but obvious contamination should be avoided.
- Make sure the interior of the restoration is clean and dry. To facilitate future removal of implant restorations, lubricate the surface of the abutment with a thin layer of water-based lubricant.
- Bleed the syringe before installation of the automix dispensing tip. Mount a new, automix dispensing tip onto the CyberTemp Long dual-barrel syringe and bleed the syringe again. This insures complete flow of base and catalyst components. Syringe a thin layer of the CyberTemp Long cement directly into the restoration. You have 45 - 60 seconds working time.
- Firmly seat the crown and either hold in place or have the patient gently close on a cotton roll until the initial set which occurs in 2 - 2.5 minutes. The initial set is characterized by a distinct gel phase. Excess cement can be removed safely after 2.5 minutes, but before final set
- Final cement set is between 4 - 5 minutes, depending on conditions.

PRECAUTION:

Though rare, some individuals may be sensitive to this product. If an allergic reaction occurs, discontinue use, remove from patient and refer patient to a physician. Important Application Information: Only a small bead of Cem-Implant around the margin is necessary to securely hold your crown. This is especially true with highly retentive restorations.

INGREDIENTS:

- Multifunctional methacrylates.
- Urethane dimethacrylate
- Polymerization activator

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