CYBERTECH POST

Instructions for use:

- Prepare the root canal by using:
 1. Gates or Largo drills to remove the gutta-percha. **(A)**
 2. Cybertech Post "Pilot" drill to enlarge the canal.
 3. Cybertech Post Cylindrical-Conical or Conical drills for preparing the orifice for housing the root pin, depending on the practitioner's preference. **(B)**
- For greater safety, use drills of progressively larger diameter in order to gradually enlarge the canal to the desired diameter.
- Make sure that you retain a sufficiently thick wall of dentine.
- Insert the metal post for taking the cast. **(C)** Send the cast with post attached to the laboratory enclosing the corresponding burn-out post.
- In case of divergent roots, enclose the corresponding Cybertech Post lock(s).
- Insert a metal post with the same diameter into the canal and make a temporary crown. Only put temporary cement around the margins of the crown and not on the post.
- Have the master cast made in the laboratory.
- Insert the burn-out post into the master cast and produce the model which is poured using conventional methods.
- Finally, try the crown for fit and cement for definitive restoration. **(D)**

Precautions to use:

- Do not use burn-out post for taking the cast. There is a difference in size between the metal and burn-out post of 2/100, the burn-out ones being smaller to allow for the thickness of the cement and to prevent excess friction between the metal and the tooth.
- Pay attention to the difference in the lengths and the diameters of the posts.
- Cybertech Post stainless steel posts must be sterilised before use.
- Cybertech Post drills and stops must be sterilised before first use and after each use.
- Cybertech Post burn-out posts and the main kits cannot be sterilised.
- A Dental dam must be used.
- Cybertech Post drills must be replaced every 40 to 60 uses.
- Cleaning and sterilisation: All the metallic parts can be cleaned with ultrasonic baths and sterilised in autoclave. Posts and instruments must be cleaned before use by dumping in a disinfecting and decontaminating solution, then sterilised (to 135 °C during 18 minutes - time of preservation) according to a validated process.





Caution!

General risks:

- This product is intended for practitioners who have received suitable training.
- If you have any doubts about using this product, consult your local supplier or the manufacturer.
- This is risk of POSTS or accessories falling into the mouth so use a Dental dam
- The choice of the size, shape and of the post material is up to the practitioner.
- The burn-out post are breakable, so never use them in the mouth.

Risks of infections : To prevent risks of infections follow recommendations concerning sterilisation.

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DE Healthcare Products
Gillingham
ME8 0SB U.K.
www.cybertechbrand.com

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