# CYBERTECH SILICONE EXTRA MATT

## Vinyl Polysiloxane Impression Material/Hydrophilic Definition


Cybertech Silicone Extra Matt Light Body & Light Body Fast is a silicone-based impression material for use in dentistry with automatic mixing device.

**Material type**  
Polyvinylsiloxane, addition-type silicone elastomer.  
light body / fast light body

– low viscosity  
Base: light green / Catalyst: white  
– ISO 4823, Type 3, low consistency

### Indications

– Wash material for the putty-wash impression technique  
– Syringe material for the simultaneous mixing technique  
– Impression material for relining

<b>Light body</b>	<b>Light body fast</b>
 	 
≤ 60 sec	≥ 120 sec
≥ 180 sec	≥ 130 sec

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### Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

### Side effects / Interactions

Polyvinylsiloxanes have a very good biological compatibility and have so far shown no harmful reactions or secondary effects on patients or dental personnel. Impression materials are intended to set in the mouth of the patient, however, the time in the mouth should not be more than double the setting time indicated. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

### Trays

Tray selection depends on impression technique and professional preference. For perfect adhesion, we recommend applying a thin layer of adhesive specified for use with polyvinylsiloxane impression materials.

### Dispensing preparation

- Remove the lid of the cartridge and discard it.
- Eject a small amount of material directly out to a pad to check that the material extrudes correctly.
- Immediately remove materials from the cartridge.
- Attach the selected mixing tip and twist a 1/4 turn clock-wise (90 °) to lock in place.
- Push the Oral Tip firmly onto the Mixing Tip – if required.
- Dispense with a smooth, steady force to mix material. Release the trigger to stop the flow.

### Mixing Tip

Do not remove the mixing tip as this acts as the cap. Store in horizontal position.. Attach a new mixing tip when the cartridge is used again ensuring that it dispenses correctly before doing so.

### Recommendations for use

-double mix technique:  
While loading the tray, the dentist may begin to apply the wash-material around the preparation. Tray loading and application of wash-material around the preparation have to be completed at the same time. Place the oral tip immersed in the material at all times to avoid trapping air. Immediately place fill tray into the patient's mouth. Press for 2-3 s and hold in position until the material is completely set.

### -two step technique-

If a 2-Step -putty /wash- technique is used, the primary impression, must be carefully cleaned and dried before use. This will help guarantee good adhesion between the wash (low viscosity material) and the material.

**For extended working time, material can be placed in the refrigerator prior to use. Check intraorally if the material is completely set before removing from the mouth.**

**Disinfection**  
The impression should be rinsed under running tap water after removal from the mouth. After rinsing, disinfection with a suitable impression material disinfectant. This will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

### Fabrication of models

The stone models can be poured after 30 min at the earliest. The impression remains dimensionally stable for at least 1 week. The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in lukewarm clear water afterwards. Cybertech impression materials are compatible with all dental stones, epoxy resin and polyurethane.

### Electro plating

Compatible with all commercially available copper and silver plating baths.

### Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent will dissolve the adhesive. These solvents should only be used in a well ventilated area. The trays can also be cleaned in an ultrasonical and sterilized if they are reusable.

### Shell life and storage

CYBERTECH SILICONE EXTRA MATT LIGHT BODY & LIGHT BODY FAST will remain usable until the expiry date provided they are stored seal and at the temperatures stated. (15-23 °C / 59-73 °F and 50 % relative air humidity) Impressions should be stored at normal room temperature, avoid exposure to heat and sun.

### Caution

Federal law restricts this device to sale by or on the order of a dentist.

### Marking

The expiry date and number are shown on the package and cartridge.

### Technical data ISO 4823-2000

Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

**MSDS** → [www.cybertechbrand.com](#)

<b>Additional products:</b>					
<b>9002833</b> Universal Autoklav Dispenser					
<b>9003270</b> Mixing canulas, yellow					
9008885 Light body fast	1		6		
9008699 Light body	1		6		

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SILICONE EXTRA MATT

UK

#### Vinyl Polysiloxane Impression Material

**Definición**
CYBERTECH SILICONE EXTRA MATT is a silicone-based impression material for use in dentistry.


**Materiale type**
Polyvinylsiloxane, addition-type silicone elastomer.

**Putty Soft**
Base: brown
Catalyst: light grey
– ISO 4823, Type 0, putty consistency

**Fast Putty Soft**
Base: brown
Catalyst: light grey
– ISO 4823, Type 0, putty consistency

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Putty Soft / Fast Putty Soft
– Primary impression in the putty-wash impression technique
– Tray material for the simultaneous mixing technique
– Tray material for the copelandhand impression technique
– Impressions for study models, orthodontic models, matrices

					
≤ 30 sec	≤ 40 sec	≥ 120 sec	≤ 30 sec	≤ 35 sec	≥ 110 sec
≥ 190 sec (3:10 min)			≥ 175 sec (2:55 min)		

**Important**

- Check intraorally if the material is completely set before removing.
- High temperatures will speed up the setting process and low temperatures will slow it down.
- For extended working time, material can be placed in the refrigerator prior to use.

**Contra-indicaciones**
Provided that the product is used as directed, there are no known contra-indications.

**Side effects / Interactions**
Polyvinylsiloxanes have very good biological compatibility and have so far shown no harmful reactions or secondary effects on patients or dental personnel. Impression materials are intended to set in the mouth of the patient, however, the time in the mouth should not be more than double the setting time indicated. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

**Trays**

Tray selection depends on impression technique and professional preference. For perfect adhesion, we recommend applying a thin layer of CYBERTECH Tray Adhesive or any other brand of adhesive specified for use with polyvinylsiloxane impression materials.

**Dosage**
Measure out equal scoops of base and catalyst.

**Attention:** Avoid interchanging the color coded measuring-scoops, each belonging to one of the components, since there is a risk of contamination!

**Mixing**

- Take the base material and the catalyst with the fingertips and start mixing. (Always wear gloves)
- Knead the material with the fingertips until the color is uniform.

**Attention:** Only use the Base and Catalyst from the same package.

**Recommendations for use**

-**double mix technique-**
While mixing, the dentist may begin to apply the wash-material around the preparation. Tray loading and application of wash-material around the preparation have to be completed at the same time. Keep the oral tip immersed in the material at all times to avoid trapping air. *Immediately place tray filled with putty body material into the patient's mouth. Press for 2-3 s and hold in position until the material is completely set.*

-**two step technique-**
If a 2-Step –putty/wash- technique is used, the primary impression made with putty material, must be carefully cleaned and dried before use. This will help guarantee good adhesion between the wash (low viscosity material) and putty tray material.

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
**Putty soft**
Basis: braun
Katalysator: hellgrau
– ISO 4823, Type 0, putty consistency

					
≤ 30 sec	≤ 40 sec	≥ 120 sec	≤ 30 sec	≤ 35 sec	≥ 110 sec
≥ 190 sec (3:10 min)			≥ 175 sec (2:55 min)		

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For extended working time, material can be placed in the refrigerator prior to use.

**Important**

Always wear gloves. Some latex gloves and surfaces contaminated by these gloves (teeth, preparations, retraction cords, etc.) may interfere with the setting process of polyvinylsiloxanes. It may be necessary to use vinyl gloves or to thoroughly wash, rinse and dry hands. To determine if your gloves will interfere with the setting process, extrude a small amount of material onto the palm of a gloved hand. Close your hand around the material and see if it sets within the time limits listed above. Compounds containing eugenol or haemostatic substances can also impede perfect setting. If (H2O2) hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with lukewarm water afterwards in order to avoid bubble formation. Intensive mixing or high temperature will speed up setting time and low temperature will slow this down.

**Check intraorally if the material is completely set before removing from the mouth.**

**Disinfection**

The impression should be rinsed under (cold) running tap water after removal from the mouth. After rinsing, disinfect with a suitable impression material disinfectant. This will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

**Manufacture of models**

Die stone models may be poured after 30 min at the earliest. The impression remains dimensionally stable for a practically unlimited period of time (at least 7 days). The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in lake-warm clear water afterwards. Cybertech impression materials are compatible with standardised dental stones, epoxy resin and polyurethane.

**Electro plating**

Compatible with commercially available copper plating and silver plating bat.

**Cleaning of trays**

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent will dissolve the Cybertech Adhesive. These solvents should only be used in a well-ventilated area. The trays can be cleaned and disinfected as usual.

**Shelf life and storage**

CYBERTECH SILICONE EXTRA MATT remains usable until the expiry date shown on the tubs provided they are sealed and stored at the temperature stated. Impressions should be stored at normal room temperature, avoid exposure to heat and sun.

**Marking**
The expiry date and number are shown on the package and cartridge.

**Technical data ISO 4823-2000**
Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

**Putty Soft**
Mixing time (15ml): 0:30 min - Time in mouth (35 °C / 95 °F): 2:00 min

**Fast Putty Soft**
Mixing time (15ml): 0:30 min - Time in mouth (35 °C / 95 °F): 1:50 min

**Caution:**

Legal Restricts this device to sale by or on the order of a dentist. Only supplied to dentists and dental laboratories or upon their instructions.

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<b>Additional products:</b>			
9008697	Putty soft		
9008698	Putty soft fast		

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Die Wahl des Löffels ist abhÄngig von der Abformtechnik oder persönlicher Präferenz (konfektionierte Löffel, „Dual-Arch“ oder individuelle Löffel). Für eine einwandfreie Haftung empfehlen wir alle Löffel mit einer dünnen Schicht © Adhesive oder mit jedem anderen, für Polyvinylsiloxane geeigneten Adhäsiv zu bestreichen.