900 8910

CYBERTECH

FILL

**DESENSITISING
MOUSSE**

Mousse Desensibilizante

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Desensibilisierungsmousse

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2x1,2ml

Syringes
Jeringas
Seringues
Spritzen
Siringhe

5

Tips
Puntas
Embouts
Spitzen
Puntali

INSTRUCTIONS FOR CYBERTECH DESENSITISING MOUSSE

Paste for remineralisation and desensitising of teeth.

The product is for use either by a dental professional in the office or provided to the patient for home treatment of dentin sensitivity.

Cybertech Desensitising Mousse provides intensive care for teeth also following a tooth whitening procedure. It remineralises the teeth and can relieve possible sensitivity. The cream has a pleasant fruity taste. Cybertech Desensitising Mousse contains Sodium fluoride, Potassium nitrate and 30 % Xylitol. The white paste can be used both in whitening trays, and directly on the teeth.

Procedure:

- Fabricate tray to meet patient's needs. (See laboratory instructions) Cybertech Desensitising Mousse may be applied in the same tray and manner as bleaching gel or into a bleaching tray if available.

Note: If large retention areas exist from prior periodontal surgery that might present problematic impression removal from the teeth, block out embrasures with Block out material

- Cybertech Desensitising Mousse is dispensed in 1.2ml unit dose syringes. A complete arch should use most of the contents of a syringe. If sensitivity is localized, place paste in affected area(s) and a small amount in the molar region to hold tray in place.

- The length and number of times a tray is worn depends on condition, patient and clinician. Recommended treatment times can range from 15 minutes to 1 hour.
- Tray should be removed at meal times. Immediately following tray removal, patient should brush and rinse teeth with water to remove all remaining paste.
- DO NOT swallow solution during rinsing. Wash tray and store in a tray storage case. Fresh paste should be used with each application.
- Paste may be warmed before placing in tray to avoid thermal sensitivity.

Laboratory Instructions:

- Pour and trim cast as you would for a standard bleaching tray. Apply 0.5-1.0mm of block-out resin to a clean, dry cast on areas in which patient is experiencing sensitivity (see Fig. 1 & 2). Also apply resin to areas where there are open retention areas, etc. For periodontal patients, clinician may determine that preoperative impressions are indicated. In such cases, reservoir spacing will be created post-op because of recontoured and apically positioned gingiva.
- Cure block out according to manufacturer's instructions in a light curing unit. Wipe off oxygen inhibition layer.
- Prepare a soft splint acc. to manufacturer's instructions.
- Trim excess tray material, allowing 4-6mm to overlap onto gingival tissue. Contour tray margin to avoid anatomical structures such as frenum attachments, bone growth, etc.

- Return tray to model; check tray extensions. Gently flame polish edges, one quadrant at a time.
- While still warm, hold periphery of each segment firmly against model for three seconds with watermoistened finger. If an area is short of the desired length, gently heat and push tray material to the desired location. If this thins the material too much, a new tray should be fabricated.

Use with Bleaching Products

Office Bleach – First Treatment

Apply Cybertech Desensitising Mousse directly to the teeth and leave for about 10 minutes, then rinse.

Home Bleach – Continuing Bleaching Treatment

After tooth whitening treatment, clean the trays and fill with cream. Wear for a minimum of 10 minutes. Tip: In case of sensitive teeth Cybertech Desensitising Mousse can be applied for 10 minutes prior to the tooth whitening procedure.

Cybertech Desensitising Mousse

Syringe with 1.2 ml paste

Ingredients: Pure water, Xylitol, Glycerine, Potassium nitrate, Cetylalcohol, Sodium hydroxide, Polyacrylic acid, Titanium dioxide, disodium hydrogen phosphate, Sodium fluoride (1450 ppm), Methyl-4-hydroxybenzoate, flavour.

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DE Healthcare Ltd